



# SINGLE-USE FLEXIBLE VIDEO BRONCHOSCOPE DISS™

DIRECT-IN-SCOPE-SUCTION





**1 DISS™ DIRECT-IN-SCOPE-SUCTION**

**2 180° ANGULATION**

**3 REDUCED OUTER DIAMETER PROVIDES IMPROVED PATIENT COMFORT**

- **STRAIGHT WORKING CHANNEL** FACILITATING EASY INSTRUMENT ACCESS
- **HD IMAGE PROCESSOR** PV300 PLUG & PLAY & RECORD STILL IMAGES & VIDEO FOOTAGE & CONNECT DIRECT TO HOSPITAL STACK SYSTEM



Model	BF220-312S Standard Deflection	BF220-31R Reverse Deflection	BF220-528S Standard Deflection	BF220-528R Reverse Deflection
Field of View	90°		90°	
Outer Diameter	4.8mm		5.2mm	
Working Channel Diameter	2.4mm		2.8mm	
Deflection Up / Down	Up / Down 180°		Up / Down 180°	
Custom Function Buttons	2		2	

Product Code	Product Description	Box Quantity
BF220-528S	Single Use Flexible Video Bronchoscope 528- Standard Deflection, OD 5.2mm, ID 2.8mm	10
BF220-312S	Single Use Flexible Video Bronchoscope 312- Standard Deflection, OD 4.8mm, ID 2.4mm	10
BF220-312R	Single Use Flexible Video Bronchoscope 312- Reverse Deflection OD 4.8mm, ID 2.4mm	10
BF220-528R	Single Use Flexible Video Bronchoscope 528- Reverse Deflection, OD 5.2mm, ID 2.8mm	10

## Customer Services

Order Direct: [office@cjmedical.com](mailto:office@cjmedical.com)

**CJ Medical**  
Charles House, Three Milestone, Truro, Cornwall, TR4 9FB  
**020 3198 2420**

@cjmedicaluk  
 [cjmedical.com](http://cjmedical.com)  
 CJ Medical UK