



# HYSTEROSCOPE SINGLE-USE

Facilitates 'See and Treat' Outpatient Hysteroscopy<sup>1</sup>



**Avoids Patient Anesthesia**




**140° Field Of View**




**HD Image**



**Optional Working Channel**



- 1 **2.9MM - 3.5MM DIAMETER - PROMOTING BEST PRACTICE NICE GUIDELINES**<sup>1</sup>
- 2 **HD IMAGE PROCESSOR** PROVIDES PREMIUM IMAGE QUALITY
- 3 **VAGINOSCOPIC APPROACH** AVOIDS CERVICAL DILATION, ANESTHESIA, AND SPECULUM USE - IMPROVED PATIENT COMFORT & OUTCOMES<sup>2</sup>
- 4 **140° FIELD OF VIEW** FACILITATES FULL VIEW OF UTERUS AND FALLOPIAN TUBES
- 5 **STRAIGHT WORKING CHANNEL** FACILITATES EASY INSTRUMENT ACCESS
-  **CONTAINS RECYCLED MATERIALS - REDUCED THEATRE EMISSIONS (KgCO2e)**

| Product Code | Product Description                                  | Box Quantity (Units) |
|--------------|------------------------------------------------------|----------------------|
| HO35         | 3.5mm Single Use Hysteroscope w/ 2mm working channel | 10                   |
| HO42         | 4.2mm Single Use Hysteroscope w/ 2mm working channel | 10                   |
| HC29         | 2.9mm Single Use Hysteroscope                        | 10                   |
| HC35         | 3.5mm Single Use Hysteroscope                        | 10                   |
| S13          | Hysteroscope Image Processor                         | 1                    |

## CUSTOMER SERVICES

ORDER DIRECT: OFFICE@CJMEDICAL.COM

**CJ Medical**

Charles House, Three Milestone, Truro, Cornwall, TR4 9FB

**020 3198 2420**

 cjmedicaluk

 cjmedical.com

 CJ Medical UK

<sup>1</sup> <https://www.nice.org.uk/guidance/gs47/chapter/Quality-statement-2-Outpatient-hysteroscopy>

<sup>2</sup> <https://www.rcog.org.uk/guidance/browse-all-guidance/green-top-guidelines/outpatient-hysteroscopy-green-top-guideline-no-59/>